

# Eukaryotic translation initiation factor 4E-binding protein 1

Catalog No: #AP80053

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP80053-1 50ug #AP80053-2 100ug #AP80053-3 1mg

## Description

Product Name	Eukaryotic translation initiation factor 4E-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	4E-BP1, Phosphorylated heat- and acid-stable protein regulated by insulin 1, PHAS-I
Accession No.	Uniprot ID:Q13541
Uniprot	Q13541
GeneID	1978;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: EIF4EBP1

Note: This product is for in vitro research use only