

14-3-3 protein sigma

Catalog No: #AP80056



Package Size: #AP80056-1 50ug #AP80056-2 100ug #AP80056-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	14-3-3 protein sigma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Stratifin,Mkrn3
Accession No.	Uniprot ID:O70456
Uniprot	O70456
GeneID	55948;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sfn

Note: This product is for in vitro research use only