

N-terminal propeptide of Collagen alpha-1(I) chain

Catalog No: #AP80058



Package Size: #AP80058-1 50ug #AP80058-2 100ug #AP80058-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	N-terminal propeptide of Collagen alpha-1(I) chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Accession No.	Uniprot ID:P02452
Uniprot	P02452
GeneID	1277;
SDS-PAGE MW	9.79kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:46-134Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PINP

Note: This product is for in vitro research use only