

## RFP-Tag Monoclonal Antibody(Mix)

Catalog No: #44063

Package Size: #44063-1 50ul #44063-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	RFP-Tag Monoclonal Antibody(Mix)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	Mix
Isotype	IgG1
Applications	WB
Specificity	The antibody detects RFP tag fusion proteins.
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	RFP-Tag
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western blotting: 1:500-1:2000

## Images



Western blot analysis of recombinant RFP protein, diluted at 1:10,000.

Note: This product is for in vitro research use only