

4F2 cell-surface antigen heavy chain

Catalog No: #AP80069



Package Size: #AP80069-1 50ug #AP80069-2 100ug #AP80069-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	4F2 cell-surface antigen heavy chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	4F2hc,Solute carrier family 3 member 2,CD98
Accession No.	Uniprot ID:Q794F9
Uniprot	Q794F9
GeneID	50567;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slc3a2

Note: This product is for in vitro research use only