

Eukaryotic translation initiation factor 4E-binding protein 1

Catalog No: #AP80073

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP80073-1 50ug #AP80073-2 100ug #AP80073-3 1mg

Description

Product Name	Eukaryotic translation initiation factor 4E-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	4E-BP1, Phosphorylated heat- and acid-stable protein regulated by insulin 1, PHAS-I
Accession No.	Uniprot ID:Q62622
Uniprot	Q62622
GeneID	116636;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Eif4ebp1

Note: This product is for in vitro research use only