

C-telopeptide of Collagen alpha-1(I) chain

Catalog No: #AP80077



Package Size: #AP80077-1 50ug #AP80077-2 100ug #AP80077-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	C-telopeptide of Collagen alpha-1(I) chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Accession No.	Uniprot ID:P02454
Uniprot	P02454
GeneID	29393;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CTXI

Note: This product is for in vitro research use only