

7-methylguanosine phosphate-specific 5'-nucleotidase

Catalog No: #AP80078



Package Size: #AP80078-1 50ug #AP80078-2 100ug #AP80078-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | 7-methylguanosine phosphate-specific 5'-nucleotidase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | 7-methylguanosine nucleotidase, Cytosolic 5'-nucleotidase 3B, Cytosolic 5'-nucleotidase III-like protein, cN-III-like protein, N(7)-methylguanylate 5'-phosphatase, Nt5c3l |
| Accession No. | Uniprot ID:Q6AYP7 |
| Uniprot | Q6AYP7 |
| GeneID | 360629; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Nt5c3b

Note: This product is for in vitro research use only