

Retinal-specific ATP-binding cassette transporter

Catalog No: #AP80081



Package Size: #AP80081-1 50ug #AP80081-2 100ug #AP80081-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Retinal-specific ATP-binding cassette transporter
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ATP-binding cassette sub-family A member 4,RIM ABC transporter,RIM protein,Stargardt disease protein,ABCR
Accession No.	Uniprot ID:P78363
Uniprot	P78363
GeneID	24;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ABCA4

Note: This product is for in vitro research use only