

Putative alpha-1-antitrypsin-related protein

Catalog No: #AP80083



Package Size: #AP80083-1 50ug #AP80083-2 100ug #AP80083-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Putative alpha-1-antitrypsin-related protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Protease inhibitor 1-like,Serpin A2,ARGS,ATR,PIL,SERPINA2P
Accession No.	Uniprot ID:P20848
Uniprot	P20848
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SERPINA2

Note: This product is for in vitro research use only