

# 7,8-dihydro-8-oxoguanine triphosphatase

Catalog No: #AP80085



Package Size: #AP80085-1 50ug #AP80085-2 100ug #AP80085-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	7,8-dihydro-8-oxoguanine triphosphatase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	2-hydroxy-dATP diphosphatase,8-oxo-dGTPase,Nucleoside diphosphate-linked moiety X motif 1,Nudix motif 1,MTH1
Accession No.	Uniprot ID:P36639
Uniprot	P36639
GeneID	4521;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:NUDT1

Note: This product is for in vitro research use only