

Aspartate aminotransferase, cytoplasmic

Catalog No: #AP80089

Package Size: #AP80089-1 50ug #AP80089-2 100ug #AP80089-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aspartate aminotransferase, cytoplasmic
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	cAspAT,Cysteine aminotransferase, cytoplasmic,Cysteine transaminase, cytoplasmic,cCAT,Glutamate oxaloacetate transaminase 1,Transaminase A
Accession No.	Uniprot ID:P00503
Uniprot	P00503
GeneID	396967;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GOT1

Note: This product is for in vitro research use only