

Aspartate aminotransferase, mitochondrial

Catalog No: #AP80090



Package Size: #AP80090-1 50ug #AP80090-2 100ug #AP80090-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aspartate aminotransferase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	mAspAT,Fatty acid-binding protein,FABP-1,Glutamate oxaloacetate transaminase 2,Kynurenine aminotransferase 4,Kynurenine aminotransferase IV,Kynurenine--oxoglutarate transaminase 4,Kynurenine--oxoglutarate transaminase IV,Plasma membrane-associated fatty acid-binding protein,FABPpm,Transaminase A
Accession No.	Uniprot ID:P12344
Uniprot	P12344
GeneID	286886;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GOT2

Note: This product is for in vitro research use only