

Alpha-2-antiplasmin

Catalog No: #AP80091



Package Size: #AP80091-1 50ug #AP80091-2 100ug #AP80091-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha-2-antiplasmin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Alpha-2-AP,Alpha-2-plasmin inhibitor,Alpha-2-PI,Serpin F2,Pli
Accession No.	Uniprot ID:Q61247
Uniprot	Q61247
GeneID	18816;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Serpinf2

Note: This product is for in vitro research use only