

ATP-binding cassette sub-family G member 2

Catalog No: #AP80097



Package Size: #AP80097-1 50ug #AP80097-2 100ug #AP80097-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | ATP-binding cassette sub-family G member 2 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | Breast cancer resistance protein,CDw338,Mitoxantrone resistance-associated protein,Placenta-specific ATP-binding cassette transporter,Urate exporter,CD338,ABCP,BCRP,BCRP1,MXR |
| Accession No. | Uniprot ID:Q9UNQ0 |
| Uniprot | Q9UNQ0 |
| GeneID | 9429; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ABCG2

Note: This product is for in vitro research use only