

5'-AMP-activated protein kinase catalytic subunit alpha-2

Catalog No: #AP80100

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP80100-1 50ug #AP80100-2 100ug #AP80100-3 1mg

Description

Product Name	5'-AMP-activated protein kinase catalytic subunit alpha-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	AMPK subunit alpha-2,Acetyl-CoA carboxylase kinase,ACACA kinase,Hydroxymethylglutaryl-CoA reductase kinase,HMGCR kinase,AMPK,AMPK2
Accession No.	Uniprot ID:P54646
Uniprot	P54646
GeneID	5563;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PRKAA2

Note: This product is for in vitro research use only