

Amyloid beta A4 protein

Catalog No: #AP80103



Package Size: #AP80103-1 50ug #AP80103-2 100ug #AP80103-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Amyloid beta A4 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	ABPP,APP,Alzheimer disease amyloid A4 protein homolog,Amyloid precursor protein,Amyloidogenic glycoprotein,AG,Beta-amyloid precursor protein
Accession No.	Uniprot ID:P08592
Uniprot	P08592
GeneID	54226;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:App

Note: This product is for in vitro research use only