

Alpha- and gamma-adaptin-binding protein p34

Catalog No: #AP80107

Package Size: #AP80107-1 50ug #AP80107-2 100ug #AP80107-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha- and gamma-adaptin-binding protein p34
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	P34
Accession No.	Uniprot ID:Q9R0Z7
Uniprot	Q9R0Z7
GeneID	171435;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Aagab

Note: This product is for in vitro research use only