

14-3-3 protein zeta/delta

Catalog No: #AP80109



Package Size: #AP80109-1 50ug #AP80109-2 100ug #AP80109-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	14-3-3 protein zeta/delta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Protein kinase C inhibitor protein 1,KCIP-1
Accession No.	Uniprot ID:P63104
Uniprot	P63104
GeneID	7534;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:YWHAZ

Note: This product is for in vitro research use only