

4F2 cell-surface antigen heavy chain

Catalog No: #AP80110



Package Size: #AP80110-1 50ug #AP80110-2 100ug #AP80110-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	4F2 cell-surface antigen heavy chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	4F2hc,4F2 heavy chain antigen,Lymphocyte activation antigen 4F2 large subunit,Solute carrier family 3 member 2,CD98,MDU1
Accession No.	Uniprot ID:P08195
Uniprot	P08195
GeneID	6520;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SLC3A2

Note: This product is for in vitro research use only