

Disintegrin and metalloproteinase domain-containing protein 17



Catalog No: #AP80114

Orders: order@signalwayantibody.com

Package Size: #AP80114-1 50ug #AP80114-2 100ug #AP80114-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Disintegrin and metalloproteinase domain-containing protein 17
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	ADAM 17, TNF-alpha convertase, TNF-alpha-converting enzyme, CD156b, TACE
Accession No.	Uniprot ID:O77636
Uniprot	O77636
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ADAM17

Note: This product is for in vitro research use only