

Putative aspartate aminotransferase, cytoplasmic 2

Catalog No: #AP80115



Package Size: #AP80115-1 50ug #AP80115-2 100ug #AP80115-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Putative aspartate aminotransferase, cytoplasmic 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Glutamate oxaloacetate transaminase 1-like protein 1, Transaminase A-like protein 1
Accession No.	Uniprot ID:Q7TSV6
Uniprot	Q7TSV6
GeneID	76615;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Got111

Note: This product is for in vitro research use only