

Alpha-2-macroglobulin

Catalog No: #AP80118

Package Size: #AP80118-1 50ug #AP80118-2 100ug #AP80118-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha-2-macroglobulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Alpha-2-M,C3 and PZP-like alpha-2-macroglobulin domain-containing protein 5,CPAMD5,FWP007
Accession No.	Uniprot ID:P01023
Uniprot	P01023
GeneID	2;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: A2M

Note: This product is for in vitro research use only