

5'-nucleotidase

Catalog No: #AP80120

Package Size: #AP80120-1 50ug #AP80120-2 100ug #AP80120-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	5'-nucleotidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	5'-NT,Ecto-5'-nucleotidase,CD73,NT5,NTE
Accession No.	Uniprot ID:Q05927
Uniprot	Q05927
GeneID	281363;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NT5E

Note: This product is for in vitro research use only