

7,8-dihydro-8-oxoguanine triphosphatase

Catalog No: #AP80121



Package Size: #AP80121-1 50ug #AP80121-2 100ug #AP80121-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	7,8-dihydro-8-oxoguanine triphosphatase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	2-hydroxy-dATP diphosphatase,8-oxo-dGTPase,Nucleoside diphosphate-linked moiety X motif 1,Nudix motif 1,Mth1
Accession No.	Uniprot ID:P53369
Uniprot	P53369
GeneID	117260;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nudt1

Note: This product is for in vitro research use only