

Aldehyde dehydrogenase family 16 member A1

Catalog No: #AP80122



Package Size: #AP80122-1 50ug #AP80122-2 100ug #AP80122-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Aldehyde dehydrogenase family 16 member A1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Accession No.	Uniprot ID:A6QR56
Uniprot	A6QR56
GeneID	506329;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ALDH16A1

Note: This product is for in vitro research use only