

ATP-binding cassette sub-family A member 7

Catalog No: #AP80123



Package Size: #AP80123-1 50ug #AP80123-2 100ug #AP80123-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ATP-binding cassette sub-family A member 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ABCA-SSN,Autoantigen SS-N,Macrophage ABC transporter
Accession No.	Uniprot ID:Q8IZY2
Uniprot	Q8IZY2
GeneID	10347;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ABCA7

Note: This product is for in vitro research use only