

Acetyl-CoA carboxylase

Catalog No: #AP80127



Package Size: #AP80127-1 50ug #AP80127-2 100ug #AP80127-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Acetyl-CoA carboxylase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	ACC
Accession No.	Uniprot ID:P11029
Uniprot	P11029
GeneID	396504;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ACAC

Note: This product is for in vitro research use only