

Neuronal acetylcholine receptor subunit alpha-7

Catalog No: #AP80129



Package Size: #AP80129-1 50ug #AP80129-2 100ug #AP80129-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neuronal acetylcholine receptor subunit alpha-7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Accession No.	Uniprot ID:P54131
Uniprot	P54131
GeneID	282178;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CHRNA7

Note: This product is for in vitro research use only