

## HSP27 Monoclonal antibody(Mix-mA)

Catalog No: #44070

Package Size: #44070-1 50ul #44070-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	HSP27 Monoclonal antibody(Mix-mA)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	Mix-mA
Isotype	IgG1
Applications	WB IHC
Species Reactivity	Hu
Specificity	Antibody can detects endogenous HSP27 protein.
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	HSP27
Other Names	HSP27
Accession No.	Swiss-Prot:P04792Gene ID:3315
Uniprot	P04792
GeneID	3315;
Calculated MW	27kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western blotting: 1:1000-1:3000

Immunohistochemistry: 1:50-1:200

## Images



Western blot analysis of Hela with HSP27 Mouse mAb diluted at 1:2,000



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using HSP27 Mouse mAb diluted at 1:200.

---

Note: This product is for in vitro research use only