

Alpha-2-macroglobulin

Catalog No: #AP80131



Package Size: #AP80131-1 50ug #AP80131-2 100ug #AP80131-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alpha-2-macroglobulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Alpha-2-M,A2m1
Accession No.	Uniprot ID:P06238
Uniprot	P06238
GeneID	100911545;24153;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:A2m

Note: This product is for in vitro research use only