

5'-AMP-activated protein kinase subunit beta-1

Catalog No: #AP80132



Package Size: #AP80132-1 50ug #AP80132-2 100ug #AP80132-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	5'-AMP-activated protein kinase subunit beta-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	AMPK subunit beta-1
Accession No.	Uniprot ID:Q5BIS9
Uniprot	Q5BIS9
GeneID	534107;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PRKAB1

Note: This product is for in vitro research use only