

3-hydroxyanthranilate 3,4-dioxygenase

Catalog No: #AP80134



Package Size: #AP80134-1 50ug #AP80134-2 100ug #AP80134-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	3-hydroxyanthranilate 3,4-dioxygenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	3-hydroxyanthranilate oxygenase,3-HAO,3-hydroxyanthranilic acid dioxygenase,HAD
Accession No.	Uniprot ID:Q0VCA8
Uniprot	Q0VCA8
GeneID	510602;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HAAO

Note: This product is for in vitro research use only