

5-hydroxytryptamine receptor 2C

Catalog No: #AP80135



Package Size: #AP80135-1 50ug #AP80135-2 100ug #AP80135-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	5-hydroxytryptamine receptor 2C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	5-HT-2C, Serotonin receptor 2C
Accession No.	Uniprot ID:Q60F97
Uniprot	Q60F97
GeneID	450240;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: HTR2C

Note: This product is for in vitro research use only