

Kynurenine/alpha-aminoadipate aminotransferase, mitochondrial

Catalog No: #AP80136

Orders: order@signalwayantibody.com

Package Size: #AP80136-1 50ug #AP80136-2 100ug #AP80136-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Kynurenine/alpha-aminoadipate aminotransferase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	KAT/AadAT,2-aminoadipate aminotransferase,2-aminoadipate transaminase,Alpha-aminoadipate aminotransferase,AadAT,Kynurenine aminotransferase II,Kynurenine--oxoglutarate aminotransferase II,Kynurenine--oxoglutarate transaminase 2,Kynurenine--oxoglutarate transaminase II
Accession No.	Uniprot ID:Q5E9N4
Uniprot	Q5E9N4
GeneID	508929;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AADAT

Note: This product is for in vitro research use only