

ATP-binding cassette sub-family G member 2

Catalog No: #AP80137



Package Size: #AP80137-1 50ug #AP80137-2 100ug #AP80137-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ATP-binding cassette sub-family G member 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Breast cancer resistance protein 1 homolog,Urate exporter,CD338,Bcrp1
Accession No.	Uniprot ID:Q80W57
Uniprot	Q80W57
GeneID	312382;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Abcg2

Note: This product is for in vitro research use only