

Adenosine receptor A2a

Catalog No: #AP80138

Package Size: #AP80138-1 50ug #AP80138-2 100ug #AP80138-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Adenosine receptor A2a
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Accession No.	Uniprot ID:P30543
Uniprot	P30543
GeneID	25369;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Adora2a

Note: This product is for in vitro research use only