

## Acyl-CoA synthetase family member 4

Catalog No: #AP80140



Package Size: #AP80140-1 50ug #AP80140-2 100ug #AP80140-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Acyl-CoA synthetase family member 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	2-aminoadipic 6-semialdehyde dehydrogenase,Protein LYS2 homolog,Putative aminoadipate-semialdehyde dehydrogenase,Acsf4,U26
Accession No.	Uniprot ID:Q80WC9
Uniprot	Q80WC9
GeneID	231326;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Aasdh

Note: This product is for in vitro research use only