

14-3-3 protein zeta

Catalog No: #AP80144



Package Size: #AP80144-1 50ug #AP80144-2 100ug #AP80144-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	14-3-3 protein zeta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	RCJMB04_11I21
Accession No.	Uniprot ID:Q5ZKC9
Uniprot	Q5ZKC9
GeneID	425619;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: YWHAZ

Note: This product is for in vitro research use only