

# Alpha-1-antiproteinase

Catalog No: #AP80148



Package Size: #AP80148-1 50ug #AP80148-2 100ug #AP80148-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Alpha-1-antiproteinase  |
| Brief Description  | Recombinant Protein   |
| Host Species       | E.coli  |
| Purification       | Greater than 90% by SDS-PAGE                                  |
| Species Reactivity | Bovine  |
| Other Names        | Alpha-1-antitrypsin,Alpha-1-proteinase inhibitor,Serpin A1,PI |
| Accession No.      | Uniprot ID:P34955   |
| Uniprot            | P34955  |
| GeneID             | 280699;   |
| Tag Info           | His   |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)        |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SERPINA1

Note: This product is for in vitro research use only