

## Alpha-1B-glycoprotein

Catalog No: #AP80149



Package Size: #AP80149-1 50ug #AP80149-2 100ug #AP80149-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Alpha-1B-glycoprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Alpha-1-B glycoprotein, C44, Liver regeneration-related protein 1
Accession No.	Uniprot ID:Q9EPH1
Uniprot	Q9EPH1
GeneID	140656;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: A1bg

Note: This product is for in vitro research use only