

Aspartate aminotransferase, mitochondrial

Catalog No: #AP80150



Package Size: #AP80150-1 50ug #AP80150-2 100ug #AP80150-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aspartate aminotransferase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	mAspAT,Fatty acid-binding protein,FABP-1,Glutamate oxaloacetate transaminase 2,Kynurenine aminotransferase 4,Kynurenine aminotransferase IV,Kynurenine--oxoglutarate transaminase 4,Kynurenine--oxoglutarate transaminase IV,Plasma membrane-associated fatty acid-binding protein,FABPpm,Transaminase A,Got-2
Accession No.	Uniprot ID:P05202
Uniprot	P05202
GeneID	14719;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Got2

Note: This product is for in vitro research use only