

Adiponectin receptor protein 1

Catalog No: #AP80152



Package Size: #AP80152-1 50ug #AP80152-2 100ug #AP80152-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Adiponectin receptor protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Progestin and adipoQ receptor family member I,PAQR1,TESBP1A,CGI-45
Accession No.	Uniprot ID:Q96A54
Uniprot	Q96A54
GeneID	51094;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ADIPOR1

Note: This product is for in vitro research use only