

Disintegrin and metalloproteinase domain-containing protein 12

Catalog No: #AP80159

Orders: order@signalwayantibody.com

Package Size: #AP80159-1 50ug #AP80159-2 100ug #AP80159-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Disintegrin and metalloproteinase domain-containing protein 12
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ADAM 12,Meltrin-alpha,MLTN,UNQ346/PRO545
Accession No.	Uniprot ID:O43184
Uniprot	O43184
GeneID	8038;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ADAM12

Note: This product is for in vitro research use only