

TNF α Monoclonal antibody(Q34)

Catalog No: #44073

Package Size: #44073-1 50ul #44073-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TNF α Monoclonal antibody(Q34)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	Q34
Isotype	IgG1
Applications	WB IHC
Species Reactivity	Hu Rt Ms
Specificity	The antibody detects endogenous TNF α protein
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	TNF α
Other Names	TNF α
Accession No.	Swiss-Prot:P01375Gene ID:7124
Uniprot	P01375
GeneID	7124;
Calculated MW	17,25kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:1000-1:3000

Immunohistochemistry: 1:50-1:200

Images



Western blot analysis of Recombinant Human TNF α Protein with TNF α Mouse mAb at 1:2,000.


mouse BV2



TNF α

Huazhong University of Science and Technology (HUST)

The picture was kindly provided by our customer



Immunohistochemical analysis of paraffin-embedded Human Tonsil Tissue using TNF α Mouse mAb diluted at 1:50

Note: This product is for in vitro research use only