

Alpha-1B-glycoprotein

Catalog No: #AP80161



Package Size: #AP80161-1 50ug #AP80161-2 100ug #AP80161-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha-1B-glycoprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Alpha-1-B glycoprotein
Accession No.	Uniprot ID:P04217
Uniprot	P04217
GeneID	1;
SDS-PAGE MW	23.65kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:281-495aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:A1BG

Note: This product is for in vitro research use only