

Kynurenine/alpha-aminoadipate aminotransferase, mitochondrial

Catalog No: #AP80162

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP80162-1 50ug #AP80162-2 100ug #AP80162-3 1mg

Description

| | |
|--------------------|---|
| Product Name | Kynurenine/alpha-aminoadipate aminotransferase, mitochondrial |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | KAT/AadAT,2-aminoadipate aminotransferase,2-aminoadipate transaminase,Alpha-aminoadipate aminotransferase,AadAT,Kynurenine aminotransferase II,Kynurenine--oxoglutarate aminotransferase II,Kynurenine--oxoglutarate transaminase 2,Kynurenine--oxoglutarate transaminase II,Kat2 |
| Accession No. | Uniprot ID:Q9WVM8 |
| Uniprot | Q9WVM8 |
| GeneID | 23923; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aadat

Note: This product is for in vitro research use only