

ATP-binding cassette sub-family B member 5

Catalog No: #AP80165

Package Size: #AP80165-1 50ug #AP80165-2 100ug #AP80165-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ATP-binding cassette sub-family B member 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ABCB5 P-gp,P-glycoprotein ABCB5
Accession No.	Uniprot ID:Q2M3G0
Uniprot	Q2M3G0
GeneID	340273;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ABCB5

Note: This product is for in vitro research use only