

Bile salt export pump

Catalog No: #AP80166



Package Size: #AP80166-1 50ug #AP80166-2 100ug #AP80166-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bile salt export pump
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	ATP-binding cassette sub-family B member 11,Sister of P-glycoprotein,Bsep,Spgp
Accession No.	Uniprot ID:O70127
Uniprot	O70127
GeneID	83569;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Abcb11

Note: This product is for in vitro research use only