

ABI gene family member 3

Catalog No: #AP80167



Package Size: #AP80167-1 50ug #AP80167-2 100ug #AP80167-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	ABI gene family member 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	New molecule including SH3,Nesh,Nesh
Accession No.	Uniprot ID:Q8BYZ1
Uniprot	Q8BYZ1
GeneID	66610;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Abi3

Note: This product is for in vitro research use only